

Human FOLH1/PSMA protein (Recombinant-Active) (N-His) (STJP012642)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells
PRODUCT PROPERTIES	
Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM MES, 150mM NaCl, 5% Trehalose, pH 5.5.
Molecular Weight	80.6 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	2346
Gene Symbol	FOLH1
UniProt ID	FOLH1_HUMAN
Immunogen Region	Lys44-Ala750
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.