

Biotinylated Human CD137/TNFRSF9/4-1BB protein (C-Avi-His) {Biotin} (STJP012644)**GENERAL INFORMATION**

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in PBS, pH 7.4.
Molecular Weight	19.9 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3604
Gene Symbol	TNFRSF9
UniProt ID	TNB9_HUMAN
Immunogen Region	Leu24-Gln186

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--