

## Biotinylated Human CD47/MER6 protein (C-Avi-His) {Biotin} (STJP012650)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

### PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 µg/ml filtered solution in 10mM Tris-Citrate, 150mM NaCl, pH 8.0.
Molecular Weight	16.4 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">961</a>
Gene Symbol	<a href="#">CD47</a>
UniProt ID	<a href="#">CD47_HUMAN</a>
Immunogen Region	Gln19-Pro139

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--