

Biotinylated Human FGL1/HFREP1 protein (C-Avi-His) {Biotin} (STJP012653)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 300mM NaCl, 10% Glycerol, 2mM EDTA, pH 8.0.
Molecular Weight	36.6 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	2267
Gene Symbol	FGL1
UniProt ID	FGL1_HUMAN
Immunogen Region	Leu23-Ile312

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--