

Human CD111/NECTIN1 protein (Recombinant-Active) (C-His) (STJP012679)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM PB, 150mM NaCl, pH 7.2.
Molecular Weight	34.99 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	5818
Gene Symbol	NECTIN1
UniProt ID	NECT1_HUMAN
Immunogen Region	Gln31-Thr334

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--