

Human Laminin 332 E8 protein (Recombinant-Active) (N-His & N-Fc & N-Fc) (STJP012687)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.2 µg/ml filtered solution in PBS, pH 7.4.
Molecular Weight	74.2 kDa & 39.8 kDa & 40.4 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	[Pro699-Gln1348] & [Gln1062-Lys1172] & [Val1074-Gln1193]
------------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--