

Biotinylated Human PCSK9 protein (C-Avi) {Biotin} (STJP012704)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 50mM HEPES, 150mM NaCl, 20% Glycerol, pH 7.4.
Molecular Weight	14 kDa & 59 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	255738
Gene Symbol	PCSK9
UniProt ID	PCSK9_HUMAN
Immunogen Region	Gln31-Gln692 (V474I, G670E)

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	----------------------------------------------------------------------------------------------------------------------