

Human PCSK9 (D374Y protein (Recombinant-Active) (C-His) (STJP012707)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells
PRODUCT PROPERTIES	
Formulation	Supplied as a 0.22 µg/ml filtered solution in 50mM HEPES, 150mM NaCl, 20% Glycerol, pH 7.4.
Molecular Weight	15-18 kDa & 60-80 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	255738
Gene Symbol	PCSK9
UniProt ID	PCSK9_HUMAN
Immunogen Region	Gln31-Gln692 (D374Y, V474I, G504R, G670E)
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.