

Human CD85k/LILRB4/ILT3 protein (Recombinant-Active) (C-His) (STJP012710)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM His-HCl, 15% Trehalose, 4% Mannitol, 0.05% Tween 80, pH 6.0.
Molecular Weight	27.1 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	11006
Gene Symbol	LILRB4
UniProt ID	LIRB4_HUMAN
Immunogen Region	Gln22-Glu259

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--