

Biotinylated Human MUC16/CA125 protein (C-Fc-Avi) {Biotin} (STJP012715)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 M µm filtered solution in 50mM Tris-HCl, 100mM Glycine, pH 7.5.
Molecular Weight	57.3 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	94025
Gene Symbol	MUC16
UniProt ID	MUC16_HUMAN
Immunogen Region	Gly12660-Met12923

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--