

Human CD316/IGSF8 protein (Recombinant-Active) (C-Fc) (STJP012726)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mg/ml filtered solution in 20mM Tris-HCl, 8% Trehalose, 50mM NaCl, 4% Mannitol, 0.02% Tween 80, pH 8.5
Molecular Weight	85.7 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	93185
Gene Symbol	IGSF8
UniProt ID	IGSF8_HUMAN
Immunogen Region	Arg28-Thr579

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--