

## Human NECTIN4/PVRL4 protein (Recombinant-Active) (C-His) (STJP012736)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

### PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM Tris-HCl, 8% Sucrose, 0.05% Tween 80, pH 8.0.
Molecular Weight	35.3 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">81607</a>
Gene Symbol	<a href="#">NECTIN4</a>
UniProt ID	<a href="#">NECT4_HUMAN</a>
Immunogen Region	Gly32-Val351

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--