

Biotinylated Human CD273/PD-L2 protein (C-His-Avi) {Biotin} (STJP012743)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells
PRODUCT PROPERTIES	
Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM PB, 150mM NaCl, pH 7.4.
Molecular Weight	25.4 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	80380
Gene Symbol	PDCD1LG2
UniProt ID	PD1L2_HUMAN
Immunogen Region	Leu20-Pro219
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.