

Human ACE2 protein (Recombinant-Active) (C-Fc) (STJP012747)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 300mM NaCl, 10% Glycerol, 100mM Glycine, pH 7.4.
Molecular Weight	110.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	59272
Gene Symbol	ACE2
UniProt ID	ACE2_HUMAN
Immunogen Region	Gln18-Ser740

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--