

## Human IL36B/IL-1F8 protein (Recombinant-Active) (Tag free) (STJP012757)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | ELISA/WB |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                     |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilized from a 0.22 Mu m filtered solution in PBS, pH 7.4.                                                                                                      |
| Molecular Weight    | 17.7 kDa                                                                                                                                                            |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt. |

### TARGET INFORMATION

|                  |             |
|------------------|-------------|
| Immunogen Region | Met1-Glu157 |
|------------------|-------------|

### ADDITIONAL INFORMATION

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|