

Human IL36A/IL-1F6 protein (Recombinant-Active) (Tag free) (STJP012766)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mg/ml filtered solution in 20mM PB, 6% Trehalose, 5mM EDTA, 0.05% Tween 80, pH 7.4.
Molecular Weight	17.1 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	27179
Gene Symbol	IL36A
UniProt ID	IL36A_HUMAN
Immunogen Region	Lys6-Phe158

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--