

Human CD257/BAFF/TNFSF13B/BLYS protein (Recombinant-Active) (N-His) (STJP012771)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 µm filtered solution in PBS, pH 7.4.
Molecular Weight	19.3 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	10673
Gene Symbol	TNFSF13B
UniProt ID	TN13B_HUMAN
Immunogen Region	Ala134-Leu285

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--