

## Human DLL3 protein (Recombinant-Active) (C-His) (STJP012780)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

### PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Histidine, 8% Sucrose, 2% Mannitol, 100mM NaCl, 0.03% Tween 80, pH 5.5.
Molecular Weight	49.4 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">10683</a>
Gene Symbol	<a href="#">DLL3</a>
UniProt ID	<a href="#">DLL3_HUMAN</a>
Immunogen Region	Ala27-Leu492

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--