

Biotinylated Mouse CD262/TRAIL-R2/DR5 protein (C-Avi-His) {Biotin} (STJP012795)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 µm filtered solution in PBS, pH 7.4.
Molecular Weight	16.4 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	21933
Gene Symbol	Tnfrsf10b
UniProt ID	TR10B_MOUSE
Immunogen Region	Asn53-Ser177

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--