

Mouse CD254/RANKL/TNFSF11 protein (Recombinant-Active) (N-His-Flag) (STJP012797)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.2 µg/ml filtered solution in 20 mM His-HCl, 100 mM NaCl, 10% Trehalose, 0.05% Tween 80, pH 5.5
Molecular Weight	23 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	Lys158-Asp316
------------------	---------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--