

Mouse TGFB1/TGF-beta-1 protein (Recombinant-Active) (Tag free) (STJP012807)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.2 Mu m filtered solution in 4 mM HCl.
Molecular Weight	12.8 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	21803
Gene Symbol	Tgfb1
UniProt ID	TGFB1_MOUSE
Immunogen Region	Ala279-Ser390

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--