

Mouse TGFB3/TGF-beta-3 protein (Recombinant-Active) (Tag free) (STJP012822)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells
PRODUCT PROPERTIES	
Formulation	Lyophilized from a 0.2 Mu m filtered solution in 50 mM Glycine-HCl, 150 mM NaCl, pH 2.5.
Molecular Weight	12.7 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene Symbol	Tgfb3 Uniprot ID= TGFB3_MOUSE "
Immunogen Region	Ala299-Ser410 (Tyr338Phe)
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.