

Mouse CD309/KDR/VEGFR-2 protein (Recombinant-Active) (C-Fc) (STJP012846)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM NaH2PO4, 150mM NaCl, 0.1M Arg, 0.1M Glu, 0.01 % Tween20, pH 7.4.
Molecular Weight	110 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	16542
Gene Symbol	Kdr
UniProt ID	VGFR2_MOUSE
Immunogen Region	Ala20-Glu762

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--