

Mouse IL2 protein (Recombinant-Active) (Tag free) (STJP012862)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.2 µm filtered solution in 20 mM Sodium Citrate, 0.2% Tween 80, pH 3.0.
Molecular Weight	17.4 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	16183
Gene Symbol	IL2
UniProt ID	IL2_MOUSE
Immunogen Region	Ala21-Gln169

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--