

Mouse NOG/Noggin protein (Recombinant-Active) (C-His) (STJP012869)

GENERAL INFORMATION

| | |
|---------------|-----------------|
| Product Type | Proteins |
| Applications | ELISA/WB |
| Host / Source | Mammalian Cells |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Formulation | Lyophilized from a 0.2 µm filtered solution in PBS, 5 mM EDTA, 5% Trehalose, pH 7.4. |
| Molecular Weight | 23.9 kDa |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

TARGET INFORMATION

| | |
|------------------|----------------------------|
| Gene ID | 18121 |
| Gene Symbol | Nog |
| UniProt ID | NOGG_MOUSE |
| Immunogen Region | Gln28-Cys232 |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|