

## Rat TGFB3/TGF-beta-3 protein (Recombinant-Active) (Tag free) (STJP012895)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

### PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.2 µm filtered solution in 50 mM Glycine-HCl, 150 mM NaCl, pH 2.5.
Molecular Weight	12.7 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">25717</a>
Gene Symbol	<a href="#">Tgfb3</a>
UniProt ID	<a href="#">TGFB3_RAT</a>
Immunogen Region	Ala301-Ser412 (Tyr340Phe)

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--