

Cynomolgus monkey CD227/MUC1 protein (Recombinant-Active) (C-Fc) (STJP012912)

GENERAL INFORMATION		
	Product Type	Proteins
	Applications	ELISA/WB
	Host / Source	Mammalian Cells
PRODUCT PROPERTIES		
	Formulation	Supplied as a 0.22 Mu m filtered solution in PBS, pH 7.4.
	Molecular Weight	40.5 kDa
	Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION		
	Immunogen Region	Leu254-Gly373
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.