

**Biotinylated Cynomolgus monkey CD137/4-1BB protein (C-His-Avi) {Biotin} (STJP012926)****GENERAL INFORMATION**

|               |                 |
|---------------|-----------------|
| Product Type  | Proteins        |
| Applications  | ELISA/WB        |
| Host / Source | Mammalian Cells |

**PRODUCT PROPERTIES**

|                     |                                                                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilized from a 0.22 Mu m filtered solution in PBS, pH 7.4.                                                                                                        |
| Molecular Weight    | 19.9 kDa                                                                                                                                                              |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

**TARGET INFORMATION**

|                  |              |
|------------------|--------------|
| Immunogen Region | Leu24-Gln186 |
|------------------|--------------|

**ADDITIONAL INFORMATION**

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|