

Biotinylated Cynomolgus monkey NECTIN4/PVRL4 protein (C-His-Avi) {Biotin} (STJP012942)**GENERAL INFORMATION**

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM PB, 100mM Arg-HCl, 10% Trehalose, 0.02% Tween 20, pH 7.5.
Molecular Weight	36.7 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	Gly32-Ser349
------------------	--------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--