

## Human CCL24/MPIF-2 protein (Recombinant) (C-His) (STJP012957)

| GENERAL INFORMATION    |                                                                                                                                                                       |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type           | Proteins                                                                                                                                                              |
| Applications           | ELISA/WB                                                                                                                                                              |
| Host / Source          | Mammalian Cells                                                                                                                                                       |
| PRODUCT PROPERTIES     |                                                                                                                                                                       |
| Formulation            | Lyophilized from a 0.22 Mg/ml filtered solution in 20mM PB, 150mM NaCl, pH 7.4.                                                                                       |
| Molecular Weight       | 11.5 kDa                                                                                                                                                              |
| Storage Instruction    | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |
| TARGET INFORMATION     |                                                                                                                                                                       |
| Gene ID                | <a href="#">6369</a>                                                                                                                                                  |
| Gene Symbol            | <a href="#">CCL24</a>                                                                                                                                                 |
| UniProt ID             | <a href="#">CCL24_HUMAN</a>                                                                                                                                           |
| Immunogen Region       | Val27-Cys119                                                                                                                                                          |
| ADDITIONAL INFORMATION |                                                                                                                                                                       |
| Note                   | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.                                                  |