

Human LEFTY2 protein (Recombinant) (N-His) (STJP012960)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mg/ml filtered solution in 20mM His-HCl, 10% Trehalose, 4% Mannitol, 100mM NaCl, 0.1% Tween 80, pH 5.5.
Molecular Weight	39.1 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	7044
Gene Symbol	LEFTY2
UniProt ID	LEFTY2_HUMAN
Immunogen Region	Phe78-Pro366

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--