

Human DLL1 protein (Recombinant) (C-Fc) (STJP012966)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mg/ml filtered solution in 20mM Tris-HCl, 50mM NaCl, 10% Trehalose, 0.05% Tween 80, 2mM EDTA, pH 8.0.
Molecular Weight	83.3 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	28514
Gene Symbol	DLL1
UniProt ID	DLL1_HUMAN
Immunogen Region	Gln18-Gly540

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--