

Human HTRA2 protein (Recombinant) (C-His) (STJP012982)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 4mM HCl.
Molecular Weight	36 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	27429
Gene Symbol	HTRA2
UniProt ID	HTRA2_HUMAN
Immunogen Region	Ala134-Glu458

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--