

## Human B4GALT3 protein (Recombinant) (C-His) (STJP012995)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

### PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 150mM NaCl, 2mM MgCl <sub>2</sub> , 10% Glycerol, pH 7.5.
Molecular Weight	41.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">8703</a>
Gene Symbol	<a href="#">B4GALT3</a>
UniProt ID	<a href="#">B4GT3_HUMAN</a>
Immunogen Region	Arg32-His393

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--