

Human B4GALT4 protein (Recombinant) (C-His) (STJP012996)

GENERAL INFORMATION		
	Product Type	Proteins
	Applications	ELISA/WB
	Host / Source	Mammalian Cells
PRODUCT PROPERTIES		
	Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 150mM NaCl, pH 7.5.
	Molecular Weight	36.77 kDa
	Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION		
	Gene ID	8702
	Gene Symbol	B4GALT4
	UniProt ID	B4GT4_HUMAN
	Immunogen Region	Gln39-Ala344
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.