

Human CTSV protein (Recombinant) (C-His) (STJP012999)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM MES, 150mM NaCl, 5% Trehalose, 5% Mannitol, 0.01% Tween 80, pH 5.5.
Molecular Weight	36.69 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	1515
Gene Symbol	CTSV
UniProt ID	CATL2_HUMAN
Immunogen Region	Val18-Val334

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--