

Native Human F8/Coagulation factor VIII protein (Native) (STJP013026)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM PB, 150mM NaCl, pH 7.2.
Molecular Weight	267 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	2157
Gene Symbol	F8
UniProt ID	FA8_HUMAN
Immunogen Region	Full length

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--