

Human F10/Factor X protein (Recombinant) (C-Fc) (STJP013032)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 M filtered solution in 20mM MES, 150mM NaCl, 0.2mM CaCl ₂ , pH 5.5.
Molecular Weight	78.2 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	2159
Gene Symbol	F10
UniProt ID	FA10_HUMAN
Immunogen Region	Asn32-Lys488

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--