

Human Latent TGF-beta-1 & GARP Complex protein (Recombinant) (N-His & C-3*Strep) (STJP013047)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in PBS, 5% Trehalose, 5% Mannitol, 0.05% Tween 80, pH 7.4.
Molecular Weight	42.1 kDa & 71.0 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	[Leu30-Ser390] & [His20-Asn627]
------------------	---------------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--