

Human IL1A/IL1F1 protein (Recombinant) (C-Fc) (STJP013062)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 µm filtered solution in 20mM Tris-HCl, 150mM NaCl, pH 7.5.
Molecular Weight	44.9 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3552
Gene Symbol	IL1A
UniProt ID	IL1A_HUMAN
Immunogen Region	Ser113-Ala271

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--