

Human APOA1 protein (Recombinant) (C-His) (STJP013069)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM PB, 15% Trehalose, 0.05% Tween 80, 50mM Methionine, pH 7.4.
Molecular Weight	30 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	335
Gene Symbol	APOA1
UniProt ID	APOA1_HUMAN
Immunogen Region	Arg19-Gln267

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--