

Human PPBP/CXCL7/NAP-2 protein (Recombinant) (C-His) (STJP013078)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM HAc-Nac, 150mM NaCl, pH 4.0.
Molecular Weight	11.3 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	5473
Gene Symbol	PPBP
UniProt ID	CXCL7_HUMAN
Immunogen Region	Ser35-Asp128

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--