

Human MMP1 protein (Recombinant) (C-His) (STJP013083)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM MES, 150mM NaCl, 2mM CaCl ₂ , 1mM DTT, 0.05%Brij35, 10% Glycerol, pH 5.0.
Molecular Weight	52.88 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	4312
Gene Symbol	MMP1
UniProt ID	MMP1_HUMAN
Immunogen Region	Phe20-Asn469

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--