

## Human CTSD/Cathepsin D protein (Recombinant) (C-His) (STJP013121)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

### PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM MES, 150mM NaCl, pH 5.5.
Molecular Weight	43.6 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">1509</a>
Gene Symbol	<a href="#">CTSD</a>
UniProt ID	<a href="#">CATD_HUMAN</a>
Immunogen Region	Leu21-Leu412

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--