

Human THBS1 protein (Recombinant) (C-10*His) (STJP013126)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 µm filtered solution in 20mM PB, 8% Trehalose, 4% Mannitol, 200mM NaCl, 0.02% Tween 80, pH 6.5.
Molecular Weight	129.2 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	7057
Gene Symbol	THBS1
UniProt ID	TSP1_HUMAN
Immunogen Region	Asn19-Pro1170

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--