

## Human GUSB protein (Recombinant) (C-His) (STJP013130)

### GENERAL INFORMATION

|               |                 |
|---------------|-----------------|
| Product Type  | Proteins        |
| Applications  | ELISA/WB        |
| Host / Source | Mammalian Cells |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Formulation         | Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 150mM NaCl, pH 8.0.  |
| Molecular Weight    | 73.6 kDa  |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

### TARGET INFORMATION

|                  |                            |
|------------------|----------------------------|
| Gene ID          | <a href="#">2990</a>       |
| Gene Symbol      | <a href="#">GUSB</a>       |
| UniProt ID       | <a href="#">BGLR_HUMAN</a> |
| Immunogen Region | Leu23-Thr651               |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|