

Human HMBS protein (Recombinant) (C-His) (STJP013133)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM PB, 150mM NaCl, 5% Trehalose, 5% mannitol, 50% Glycerol, 0.1% Tween 80, pH 7.4.
Molecular Weight	40.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3145
Gene Symbol	HMBS
UniProt ID	HEM3_HUMAN
Immunogen Region	Ser2-His361

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--