

Human FBP1 protein (Recombinant) (C-His) (STJP013149)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 200mM NaCl, 1mM DTT, 1mMEDTA, 10% Glycerol, pH 8.0.
Molecular Weight	37.8 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	2203
Gene Symbol	FBP1
UniProt ID	F16P1_HUMAN
Immunogen Region	Ala2-Gln338

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--