

## Human CD106/VCAM1 protein (Recombinant) (C-Fc) (STJP013231)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

### PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 µm filtered solution in PBS, pH 7.4.
Molecular Weight	101.2 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">7412</a>
Gene Symbol	<a href="#">VCAM1</a>
UniProt ID	<a href="#">VCAM1_HUMAN</a>
Immunogen Region	Phe25-Glu698

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	----------------------------------------------------------------------------------------------------------------------