

Human DTYMK protein (Recombinant) (C-His) (STJP013259)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM PB, 150mM NaCl, 10% Sucrose, 0.05% Tween 80, 20% Glycerol, 1mM EDTA, pH 8.0.
Molecular Weight	24.9 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	1841
Gene Symbol	DTYMK
UniProt ID	KTHY_HUMAN
Immunogen Region	Met1-Lys212

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--